



## Ellie Recliner

Model: ELI12A00R

31.5 W 33.75 D 46 H 23 SW 19.25 SD 20 SH 26.5 AH

## **Features**



**Durable Finish** 

Protea™ doesn't have any of the drawbacks of wood. The frames are moisture impervious, graffiti-resistant, easy to clean and maintenance-free.



Trendelenburg Mechanism

Engaged by staff only, the mechanism activates from either side when the chair is fully reclined to 72" and pressure is applied to the back.



Optional Upholstered Arm Panels

The Upholstered Arm Panels are available and optional.



Central Locking System

Locks all four casters with a single control.



Optional Transfer Arm

An optional Transfer Arm, available on the arm opposite the Ellie Recliner's lever, can be raised and lowered for easy patient transfers.



Optional Adjustable Headrest

An adjustable headrest is available and is optional.



**Inclined Footrest** 

The footrest tucks under the seat of the Ellie Recliner for better positioning from sit to stand.



**Push Bar** 

The push bar aids staff in positioning the recliner, and is made with Protea.



Replaceable Components

The Ellie back, seat, continuous footrest, arm panels and caster with bumper are all replaceable components.



**Safety Footrest** 

Safety footrest collapses upon overload.



Optional IV Pole Holder

An optional IV Pole Holder is available and must be attached on the side opposite the Ellie Recliner's lever.



Optional Foley Bag Hook

An optional Foley Bag Hook is available and must be added to the side opposite the Ellie Recliner's lever.



Optional Heat/ Massage

Heat and Massage is an optional feature, only available if Upholstered Arm Panels are used.



Optional Tablet Table

Available as an add-on. Not compatible with the Transfer Arm option.



Optional Side Folding Table

Available as an add-on. Not compatible with the Transfer Arm option.



6160 Peachtree Dunwoody Rd. Building C, Atlanta, GA 30328

1-877-MY-KWALU **kwalu.com**